

Expect More.


code[®]

CERNER Integration Guide

CR1100 & CR1500 & CR2700

- Not Compatible with Other Models -

VERSION 03 -RELEASE DATE MARCH 2020



codecorp.com



Configuration Guide




[codecorporation](https://www.youtube.com/codecorporation)

Prepare your reader by scanning the following barcode



CC004489_01

Program your reader by scanning the following barcode



CC004463_v1926

FOR DEFAULT CONFIGURATIONS - STOP HERE

Cerner Custom Settings

Scan the appropriate barcode to enable the desired software specific formatting:

CareAware/Cerner HS/Bridge
(ALT + Scanned Data + TAB)
(Default)



CC004446_01

Farpoint/Pathnet/Blood Bank
(Scanned Data + Enter)



CC004447_01

Scan the appropriate barcode to toggle lowercase Codabar:

Lowercare Codabar Disabled
(Default)



CC004449_01

Lowercare Codabar Enabled



CC004448_01

Scan the appropriate barcode to toggle ISBT128 Concatenation:

ISBT 128 Concatenation Disabled
(Default)



CC004451_01

ISBT 128 Concatenation Enabled



CC004450_01

Scan the appropriate barcode to toggle Pharmacy conversion of <GS> to "I":

Pharmacy Character Conversion Disabled
(Default)



CC004453_01

Pharmacy Character Conversion Enabled



CC004452_01